

Model : FE303XP



* controller panel is optional

Custom multiple output power supply for "TFE-type scanning electron microscope" .(High performance type)

Application

- Scanning electron microscope
- Electron probe micro analyzer
- FE-Auger electron microscopy

Product specifications

Input voltage

90V~240VAC single phase 1.5A(50,60HZ)

Accelerator supply (Referenced to GND)

Output voltage: -1kV~-30kV
 Absolute voltage accuracy: Less than 0.1%
 Set ability: 16bit (0.5V LSB)
 Output current: 200 μ A max
 Ripple noise: Less than 2ppm 60mVp-p@30kV
 Stability: 15ppm/1hr after 1hour warm up
 temperature coefficient: 10ppm/ $^{\circ}$ C

Filament supply (Referenced to Accelerator)

Output voltage: 0~5V
 Set ability: 16bit (1mA LSB)
 Output current: 0~3.2A
 Ripple noise: Less than 5mA p-p @3A
 Absolute current accuracy: \pm 0.01A
 Stability: 20ppm/1hour @ 3.0A
 temperature coefficient: 15ppm/ $^{\circ}$ C

Suppressor supply (Referenced to Accelerator)

Output voltage: -0.03kV~ -0.6kV
 Absolute voltage accuracy: Less than 1.0% @300V
 Set ability: 14bit (0.11V LSB)
 Output current: 150 μ A max
 Ripple noise: Less than 20mVp-p
 Stability: 50ppm/1hr after 1hour warm up
 temperature coefficient: 25ppm/ $^{\circ}$ C

Extractor supply (Referenced to Accelerator)

Output voltage: +0.4kV~+8.0kV
 Absolute voltage accuracy: Less than 0.1%
 Set ability: 16bit (0.16V LSB)
 Output current: 200 μ A max
 Ripple noise: Less than 30mVp-p @4.0kV
 Stability: 15ppm/1hr after 1hour warm up
 temperature coefficient: 10ppm/ $^{\circ}$ C

Electric static lens supply (Referenced to Accelerator)

Output voltage: -0.3K~-15kV
 Absolute voltage accuracy: Less than 0.1%
 Set ability: 16bit (0.1V LSB)
 Output current: 180 μ A max
 Ripple noise: Less than 100mVp-p
 Stability: 20ppm/1hr after 1hour warm up
 temperature coefficient: 10ppm/ $^{\circ}$ C

optional

The other specifications

All the outputs provide with the over voltage protection and the over current protection.

External control: Optical isolated RS232C

Interlock: vacuum, thermo, HV-connection

HV connector: Customer specification

Output monitor: Accelerating voltage, Lens voltage

Storage temp range: -5 $^{\circ}$ C~ +40 $^{\circ}$ C

Operating temperature

Limit: +10 $^{\circ}$ C~ +40 $^{\circ}$ C

Humidity: 80% or less

Insulation method: Air insulation (one molding)

Externals size: 480mm(W) \times 450(D) \times 249 (H)

Weight: 33kg

The model "FE303XP" is an integrated multiple output high voltage power supply specifically developed for TFE type scanning electron microscope.

* Please ask details.
 e-mail info@futex.jp

FUTEX CO., LTD. <http://www.futex.jp>
 TEL:+81-42-549-2888 FAX:+81-42-549-2900
 2-28-3 Fukujima-cho Akishima Tokyo, Japan

Update: Apr,27, 2012